

Product Datasheet

Keratin, type II cytoskeletal 4 antibody (FITC) (orb241681)

Description	Keratin, type II cytoskeletal 4 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Keratin, type II cytoskeletal 4 protein (1-534AA)
Target	KRT4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19013
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.