



Product Datasheet

TBX1 antibody (orb2417)


Description	TBX1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human TBX1 (165-270/398 aa)
Target	TBX1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	O43435
Dilution Range	WB=1:500-2000, Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.



Sample:
Muscle
(Mouse)
Lysate at 40
ug P...



Sample: Brain
(Rat) Lysate
at 40 ug
Prim...



Sample:
Muscle
(Mouse)
Lysate at 40
ug P...