

Product Datasheet

Na(+)/H(+) exchange regulatory cofactor NHE-RF3 antibody (Biotin) (orb241704)

Description	Na(+)/H(+) exchange regulatory cofactor NHE-RF3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Na(+)/H(+) exchange regulatory cofactor NHE-RF3 protein (1-519AA)
Target	PDZK1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5T2W1
Expiration Date	12 months from date of receipt.