

Product Datasheet

N-acetylgalactosamine-6-sulfatase antibody (Biotin) (orb241708)

Description	N-acetylgalactosamine-6-sulfatase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human N-acetylgalactosamine-6-sulfatase protein (27-522AA)
Target	GALNS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P34059
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.