



Product Datasheet

PTH antibody (orb241719)


Description	PTH antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Parathyroid hormone protein (32-115AA)
Target	PTH
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01270
Dilution Range	IHC-P:1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-3(Lane 1), C...



Immunohistochemical staining of human pa...



Immunohistochemical staining of human pr...