



Product Datasheet

Protein-lysine 6-oxidase antibody (orb241751)


Description	Protein-lysine 6-oxidase antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Protein-lysine 6-oxidase protein (174-417AA)
Target	LOX
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P28300
Dilution Range	WB:1:500-1:5000, IHC-P:1:200-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse li...



Immunohistochemical staining of human ki...



Western blot analysis of A549 whole cell...