

Product Datasheet

Rho guanine nucleotide exchange factor 18 antibody (Biotin) (orb241768)

Description	Rho guanine nucleotide exchange factor 18 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Rho guanine nucleotide exchange factor 18 protein (159-551AA)
Target	ARHGEF18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6ZSZ5
Expiration Date	12 months from date of receipt.