

## Product Datasheet

### Calcitonin Antibody (orb24179)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | RABBIT ANTI CALCITONIN          |
| <b>Species/Host</b>        | Rabbit                          |
| <b>Tested Applications</b> | RIA                             |
| <b>Target</b>              | Calcitonin                      |
| <b>Preservatives</b>       | No Preservative                 |
| <b>Form/Appearance</b>     | Lyophilized                     |
| <b>Note</b>                | For research use only           |
| <b>Clonality</b>           | Polyclonal                      |
| <b>Purity</b>              | Neat                            |
| <b>Source</b>              | Rabbit                          |
| <b>Expiration Date</b>     | 12 months from date of receipt. |