

Product Datasheet

threonine-protein phosphatase 2A catalytic subunit beta isoform antibody (orb241791)

Description	threonine-protein phosphatase 2A catalytic subunit beta isoform antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A catalytic subunit beta isoform protein (1-309AA)
Target	PPP2CB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62714
Dilution Range	IHC-P:1:100-500, IF:1:50-200
Expiration Date	12 months from date of receipt.