

## Product Datasheet

# threonine-protein phosphatase 2A catalytic subunit beta isoform antibody (Biotin) (orb241792)

|                            |                                                                                                            |
|----------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | threonine-protein phosphatase 2A catalytic subunit beta isoform antibody (Biotin)                          |
| <b>Species/Host</b>        | Rabbit                                                                                                     |
| <b>Reactivity</b>          | Human                                                                                                      |
| <b>Conjugation</b>         | Biotin                                                                                                     |
| <b>Tested Applications</b> | ELISA                                                                                                      |
| <b>Immunogen</b>           | Recombinant Human Serine/threonine-protein phosphatase 2A catalytic subunit beta isoform protein (1-309AA) |
| <b>Target</b>              | PPP2CB                                                                                                     |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol                                             |
| <b>Form/Appearance</b>     | Liquid                                                                                                     |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                              |
| <b>Note</b>                | For research use only                                                                                      |
| <b>Isotype</b>             | IgG                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                 |
| <b>Uniprot ID</b>          | <b>P62714</b>                                                                                              |
| <b>Dilution Range</b>      | ELISA                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                            |