

Product Datasheet

Stress-induced-phosphoprotein 1 antibody (orb241804)

Description Stress-induced-phosphoprotein 1 antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen Recombinant Human Stress-induced-phosphoprotein 1 protein (1-543AA)

Target STIP1

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID **P31948**

Dilution Range WB:1:500-1:5000, IHC-P:1:20-1:200


Expiration Date 12 months from date of receipt.



Western blot analysis of Hela whole cell...



Immunohistochemical staining of human sm...



Immunohistochemical staining of human ov...