

Product Datasheet

Tax1-binding protein 3 antibody (HRP) (orb241815)

Description	Tax1-binding protein 3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Tax1-binding protein 3 protein (1-124AA)
Target	TAX1BP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14907
Expiration Date	12 months from date of receipt.