



## Product Datasheet

# Type I inositol 1, 4, 5-trisphosphate 5-phosphatase antibody (orb241872)


|                            |                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------|
| <b>Description</b>         | Type I inositol 1, 4, 5-trisphosphate 5-phosphatase antibody                          |
| <b>Species/Host</b>        | Rabbit                                                                                |
| <b>Reactivity</b>          | Human                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                          |
| <b>Tested Applications</b> | ELISA, IF, IHC                                                                        |
| <b>Immunogen</b>           | Recombinant Human Type I inositol 1,4,5-trisphosphate 5-phosphatase protein (1-409AA) |
| <b>Target</b>              | INPP5A                                                                                |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol                        |
| <b>Form/Appearance</b>     | Liquid                                                                                |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C.<br>Avoid repeated freeze.                      |
| <b>Note</b>                | For research use only                                                                 |
| <b>Isotype</b>             | IgG                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                            |
| <b>Uniprot ID</b>          | <b>Q14642</b>                                                                         |
| <b>Dilution Range</b>      | IHC-P:1:20-1:200, IF:1:50-1:200                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                       |



Immunofluorescence analysis of HepG2 cel...



4



5-trisphosphate 5-phosphatase antibody