

Product Datasheet

Alkaline phosphatase, tissue-nonspecific isozyme antibody (Biotin) (orb241877)

Description	Alkaline phosphatase, tissue-nonspecific isozyme antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Alkaline phosphatase, tissue-nonspecific isozyme protein (20-503AA)
Target	Alpl
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09242
Expiration Date	12 months from date of receipt.