

Product Datasheet

Forkhead box protein P3 antibody (FITC) (orb241885)

Description	Forkhead box protein P3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Forkhead box protein P3 protein (1-260AA)
Target	FOXP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BZS1
Expiration Date	12 months from date of receipt.