



Product Datasheet

Bile acid R antibody (orb241905)


Description	Bile acid R antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Bile acid receptor protein (1-476AA)
Target	NR1H4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96RI1
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse stomach t...



Immunohistochemical staining of human li...



Western blot analysis of HEK293 whole ce...