

Product Datasheet

Dr1-associated corepressor antibody (orb241909)

Description	Dr1-associated corepressor antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Dr1-associated corepressor protein (4-198AA)
Target	DRAP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14919
Expiration Date	12 months from date of receipt.