



Product Datasheet

Estrogen R antibody (orb241921)


Description	Estrogen R antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IHC, WB
Immunogen	Recombinant Human Estrogen receptor protein (10-591AA)
Target	ESR1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03372
Dilution Range	WB:1:1000-1:5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
MCF-7 (1.1...



Western blot aanalysis
of Hela whole cel...



Immunohistochemistry
of paraffin-embedde...