


Product Datasheet

Heterogeneous nuclear ribonucleoprotein H antibody (orb241954)

Description	Heterogeneous nuclear ribonucleoprotein H antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Heterogeneous nuclear ribonucleoprotein H protein (2-216AA)
Target	HNRNPH1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31943
Expiration Date	12 months from date of receipt.



Western blot analysis of EC109(Lane 1), ...