

Product Datasheet

Heterogeneous nuclear ribonucleoproteins A2/B1 antibody (HRP) (orb241961)

Description	Heterogeneous nuclear ribonucleoproteins A2/B1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Heterogeneous nuclear ribonucleoproteins A2/B1 protein (1-249AA)
Target	HNRNPA2B1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22626
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.