



Product Datasheet

Muscleblind-like protein 1 antibody (orb241978)

Description	Muscleblind-like protein 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Muscleblind-like protein 1 protein (1-382AA)
Target	MBNL1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NR56
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of paraffin-embedde...



Immunohistochemistry
of paraffin-embedde...