

Product Datasheet

Nuclear R subfamily 1 group I member 3 antibody (orb241986)

Description Nuclear R subfamily 1 group I member 3 antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IHC

Immunogen Recombinant Human Nuclear receptor subfamily 1 group I member 3 protein (2-351AA)

Target NR1I3

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C.
Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID [Q14994](#)

Dilution Range IHC-P:1:20-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...