



Product Datasheet

DHRS7 antibody (orb2420)

Description	DHRS7 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human DHRS7 (61-160/339 aa)
Target	DHRS7
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35 kDa
Uniprot ID	Q9Y394
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse li...



Immunohistochemical staining of rat skin...