

Product Datasheet

Small nuclear ribonucleoprotein Sm D2 antibody (FITC) (orb242035)

Description	Small nuclear ribonucleoprotein Sm D2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Small nuclear ribonucleoprotein Sm D2 protein (1-118AA)
Target	SNRPD2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62316
Expiration Date	12 months from date of receipt.