

Product Datasheet

Transcription intermediary factor 1-beta antibody (FITC) (orb242063)

Description	Transcription intermediary factor 1-beta antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Transcription intermediary factor 1-beta protein (22-291AA)
Target	TRIM28
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13263
Expiration Date	12 months from date of receipt.