



Product Datasheet

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 1 antibody [Out of stock] (orb242109)

Description	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 1 antibody [Out of stock]
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 1 protein (1-70AA)
Target	NDUFA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15239
Dilution Range	IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of HeLa cells...



Immunohistochemistry of paraffin-embedde...