

Product Datasheet

Corticosteroid-binding globulin antibody (FITC) (orb242119)

Description	Corticosteroid-binding globulin antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Corticosteroid-binding globulin protein (23-405AA)
Target	SERPINA6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08185
Expiration Date	12 months from date of receipt.