



Product Datasheet

ATP-sensitive inward rectifier potassium channel 1 antibody (orb242129)

Description	ATP-sensitive inward rectifier potassium channel 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human ATP-sensitive inward rectifier potassium channel 1 protein (178-391AA)
Target	KCNJ1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P48048
Dilution Range	WB:1:1000-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2(Lane 1), ...



Immunofluorescence analysis of HeLa cell...