



## Product Datasheet

### FKBP6 antibody (orb24214)


|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | FKBP6 antibody                                                                           |
| <b>Species/Host</b>        | Rabbit                                                                                   |
| <b>Reactivity</b>          | Human                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | ELISA, IF, IHC, WB                                                                       |
| <b>Immunogen</b>           | Recombinant Human Inactive peptidyl-prolyl cis-trans isomerase FKBP6 protein (101-327AA) |
| <b>Target</b>              | FKBP6                                                                                    |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                              |
| <b>Form/Appearance</b>     | Liquid                                                                                   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                            |
| <b>Note</b>                | For research use only                                                                    |
| <b>Isotype</b>             | IgG                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                               |
| <b>Uniprot ID</b>          | <b>O75344</b>                                                                            |
| <b>Dilution Range</b>      | IHC-P:1:20-200, IF:1:50-200                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |



Immunofluorescence analysis of MCF-7 cel...



Immunohistochemical staining of human co...



Immunohistochemical staining of human te...