

Product Datasheet

FKBP6 antibody (HRP) (orb24215)

Description	FKBP6 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Inactive peptidyl-prolyl cis-trans isomerase FKBP6 protein (101-327AA)
Target	FKBP6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O75344
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.