

Product Datasheet

AP-2 complex subunit mu antibody (FITC) (orb242151)

Description	AP-2 complex subunit mu antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human AP-2 complex subunit mu protein (1-435AA)
Target	AP2M1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96CW1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.