



Product Datasheet

Aquaporin 1 antibody (orb242161)


Description	Aquaporin 1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	Recombinant Human Aquaporin-1 protein (220-269AA)
Target	AQP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P29972
Dilution Range	WB:1:1000-1:5000, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of HepG2 cell...



Western blot analysis of Rat heart tiss...



Immunoprecipitation analysis of Mouse sk...