

Product Datasheet

ATP-dependent RNA helicase DDX19A antibody (Biotin) (orb242166)

Description ATP-dependent RNA helicase DDX19A antibody (Biotin)

Species/Host Rabbit

Reactivity Human

Conjugation Biotin

Tested Applications ELISA

Immunogen Recombinant Human ATP-dependent RNA helicase DDX19A protein (1-478AA)

Target DDX19A

Preservatives Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.


Note For research use only

Isotype IgG

Clonality Polyclonal

Uniprot ID **Q9NUU7**

Expiration Date 12 months from date of receipt.



Western blot analysis of mouse gonad tis...