

Product Datasheet

Cation-dependent mannose-6-phosphate R antibody (FITC) (orb242171)

Description	Cation-dependent mannose-6-phosphate R antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Cation-dependent mannose-6-phosphate receptor protein (1-277AA)
Target	M6PR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20645
Expiration Date	12 months from date of receipt.