

Product Datasheet

Cytochrome b antibody (HRP) (orb242179)

Description	Cytochrome b antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Cytochrome b5 protein (2-134AA)
Target	CYB5A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00167
Expiration Date	12 months from date of receipt.