

Product Datasheet

Fatty acid-binding protein antibody (Biotin) (orb242198)

Description	Fatty acid-binding protein antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Fatty acid-binding protein, heart protein (2-133AA)
Target	FABP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05413
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.