


Product Datasheet

Fatty acid-binding protein antibody (HRP) (orb242205)

Description	Fatty acid-binding protein antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Fatty acid-binding protein, liver protein (1-127AA)
Target	FABP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07148
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): recom...