

Product Datasheet

Gamma-aminobutyric acid R-associated protein-like 2 antibody (orb242213)

Description Gamma-aminobutyric acid R-associated protein-like 2 antibody

Species/Host Rabbit

Reactivity Human, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen Recombinant Human Gamma-aminobutyric acid receptor-associated protein-like 2 protein (1-117AA)

Target GABARAPL2

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C.
Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID **P60520**

Dilution Range WB:1:500-1:5000, IHC-P:1:1000-1:2000


Expiration Date 12 months from date of receipt.



Western blot analysis of 293T whole cell...



Immunohistochemical staining of human heart...



Immunohistochemical staining of human prostate...