



Product Datasheet

Mitochondrial brown fat uncoupling protein 1 antibody (orb242229)


Description	Mitochondrial brown fat uncoupling protein 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Mitochondrial brown fat uncoupling protein 1 protein (2-307AA)
Target	UCP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25874
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ad...



Immunohistochemical staining of human br...



Immunohistochemical staining of human co...