

Product Datasheet

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 2 antibody (HRP) [Out of stock] (orb242252)

Description	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 2 antibody (HRP) [Out of stock]
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 2 protein (4-99AA)
Target	NDUFA2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43678
Expiration Date	12 months from date of receipt.