

Product Datasheet

FMO4 antibody (FITC) (orb24226)

Description	FMO4 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Dimethylaniline monooxygenase [N-oxide-forming] 4 protein (1-300AA)
Target	FMO4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31512
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.