

Product Datasheet

NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 7 antibody (Biotin) [Out of stock] (orb242262)

Description	NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 7 antibody (Biotin) [Out of stock]
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 7 protein (2-137AA)
Target	NDUFB7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P17568
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.