



Product Datasheet

Nuclear pore glycoprotein 210 antibody (orb242281)


Description	Nuclear pore glycoprotein 210 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Nuclear pore membrane glycoprotein 210 protein (28-238AA)
Target	NUP210
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TEM1
Dilution Range	IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of paraffin-embedde...



Immunofluorescent
analysis of MCF-7
cell...



Immunohistochemistry
of paraffin-embedde...