

Product Datasheet

Secretory carrier-associated protein 3 antibody (FITC) (orb242307)

Description	Secretory carrier-associated protein 3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Secretory carrier-associated membrane protein 3 protein (2-170AA)
Target	SCAMP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14828
Expiration Date	12 months from date of receipt.