



Product Datasheet

FRS3 antibody (orb24232)


Description	FRS3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Fibroblast growth factor receptor substrate 3 protein (101-400AA)
Target	FRS3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	043559
Dilution Range	IHC-P:1:20-200, IF:1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ce...



Immunofluorescence analysis of U251 cell...



Immunohistochemical staining of human sk...