

Product Datasheet

V-type proton ATPase subunit E 1 antibody (orb242342)

Description	V-type proton ATPase subunit E 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human V-type proton ATPase subunit E 1 protein (2-224AA)
Target	ATP6V1E1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P36543
Dilution Range	WB:1:1000-1:5000, IF:1:50-1:200
Expiration Date	12 months from date of receipt.

