

Product Datasheet

Vascular cell adhesion protein 1 antibody (orb242379)

Description	Vascular cell adhesion protein 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Vascular cell adhesion protein 1 protein (25-698AA)
Target	VCAM1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19320
Expiration Date	12 months from date of receipt.