


Product Datasheet

U7 snRNA-associated Sm-like protein LSM10 antibody (orb242395)

Description	U7 snRNA-associated Sm-like protein LSM10 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human U7 snRNA-associated Sm-like protein LSM10 protein (1-123AA)
Target	LSM10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q969L4
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.