



Product Datasheet

DHRS4L2 antibody (orb2424)

Description	DHRS4L2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human DHRS4L2 (21-120/230 aa)
Target	DHRS4L2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	Q6PKH6
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Rat kidn...



Immunohistochemical staining of colon ca...