

Product Datasheet

Sulfolobus islandicus Chromatin protein Cren7 antibody (orb242402)

Description	Sulfolobus islandicus Chromatin protein Cren7 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Sulfolobus islandicus Chromatin protein Cren7 protein (1-60AA)
Target	creN7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	C3MPN0
Expiration Date	12 months from date of receipt.